



CHAIRSIDE®
Block Out Material
IDEAL FOR BLOCKING OUT UNDERCUTS
WHILE PROCESSING OVERDENTURES CHAIRSIDE

4 x 1.2ml syringes **\$48** #9574

The latest advancement in the CHAIRSIDE Product Portfolio; - Introducing CHAIRSIDE Block Out Material from ZEST Anchors, a light-cured protective material that removes undercuts while processing overdentures chair side.

CHAIRSIDE Denture Prep and Polish Kit

- A comprehensive kit consisting of Recess, Trim, Undercut, Grind and Vent Burs, as well as a Polisher, all of which are designed to address the most frequent overdenture preparation requirements.
- A unique to the market CHAIRSIDE Recess Bur, specifically designed by ZEST Anchors, that quickly and easily prepares the exact size recess for the LOCATOR and SATURNO Denture Attachment Housings (Denture Caps).




#9582
CHAIRSIDE Undercut Bur, Vent Bur, Trim Bur, Grind Bur, Recess Bur and Polishing Bur

\$320

CHAIRSIDE®
Denture Removal Tool

#9549
CHAIRSIDE Denture Removal Tool

\$58 (10/pcs.)



GIVE YOUR PATIENTS THE GIFT OF EASY DENTURE REMOVAL



CHAIRSIDE Attachment Processing Material

- Cost Efficient
- No primer required for reduced chair time
- Shorter mixing tip reduces material waste and cost per case*
- Simple Application
- Ideal viscosity for accurate application
- Bonds to itself for efficient void fill and touch-ups
- Dual cure capability provides convenient self or light cure flexibility
- Patient Friendly
- Low curing temperature ensures patient comfort
- Odorless and tasteless to increase patient acceptance



\$312

#9565

CHAIRSIDE® Attachment Processing Material
(18ml cartridge and 10 - Mixing Tips)



\$105

#9566-P

CHAIRSIDE® Attachment Processing Material 8g Syringe
(15 mixing and angled tips included in each package)

RELINING CHAIR SIDE? TRY CHAIRSIDE® SOFT RELINE MATERIAL

Ideal for immediate dentures following extractions, surgery or implant procedures. Silicone based material does not increase temperature while setting, is odorless, tasteless and does not harden over time. Long lasting smooth finish prevents contamination and odor. Stain-resistant, high tear strength and strong adhesion for long lasting durability up to 18 months.



CHAIRSIDE®
Soft Reline Material

\$250

#9573

