

Original Legacy Locator Components

- 8510 Locator Denture Cap w/proc male (4)
- 8510-10 Locator Denture Cap w/proc male (10)
- 8515 Locator Black Processing Male (4)
- 8515-20 Locator Black Processing Male (20)

- 8524 Locator Replace Male Clear- Increased (4)
- 8524-20 Locator Replace Male Clear- Increased (20)
- 8527 Locator Replace Male Pink- Medium (4)
- 8527-20 Locator Replace Male Pink- Medium(20)
- 8529 Locator Replace Male Blue- Decreased (4)
- 8529-20 Locator Replace Male Blue- Decreased (20)

- 8547 Locator Replace Male Green- Increased (4)
- 8547-20 Locator Replace Male Green- Increased (20)
- 8915 Locator Replace Male Orange- Medium (4)
- 8915-20 Locator Replace Male Orange- Medium (20)
- 8548 Locator Replace Male Red- Decreased (4)
- 8548-20 Locator Replace Male Red-Decreased (20)

- 8558 Locator Replace Male Gray- Zero Ret. (4)
- 8558-20 Locator Replace Male Gray- Zero Ret. (20)
- 8559 Locator Replace Male Purple-Strong (4)
- 8026 Locator Yellow Bar Processing Males (4)
- 8026-20 Locator Yellow Bar Processing Males (20)

- 8587-2 Locator Bar Female (2) (Bar) threaded 2.56mm
- 8587-10 Locator Bar Female (10) (Bar) threaded 2.56mm

- 8586-2 Locator Cast-to Female (2) (Bar)
- 8586-10 Locator Cast-to Female (10) (Bar)
- 8588-2 Locator Laser TI Bar Female (2)
- 8588-10 Locator Laser TI Bar Female (10)
- 8589-2 Locator Bar Female (2) threaded 2mm
- 8589-10 Locator Bar Female (10) threaded 2mm
- 8590-2 Locator Laser SS Bar Female (2)
- 8590-10 Locator Laser SS Bar Female (10)

- 8014 Locator 2mm Castable Threaded Inserts (10)
- 8016 Locator Bar Drill Tap Holder
- 9102 Locator 1.7mm Drill (2.0mm)
- 9104 Locator Bar Tap (2.0mm)
- 8922 Locator Spot Face Diamond Bur
- 8922 Locator Countersink Diamond Bur
- 8924 Locator Pilot Drill
- 8926 Locator Square Driver (15mm)
- 8927 Locator Square Driver (21mm)
- 8913 Locator Insert Driver (Latch 23mm)
- 8914 Locator Insert Driver (Latch 29mm)
- 8260 Locator 35n Square Driver (15mm)
- 8280 Locator 35n Square Driver (21mm)

- 1747 Locator ERA Retrofit Kit (2)
- 8703 Locator Implant Abutments (ALL)*
- 8909-2 Locator for Multi-Unit Abutment TI (2)

- 9566-P Chairside processing material (8 g)
- 9582 Chairside denture prep and polish kit
- 9120 Overdenture implant education model



8510 Denture Cap w/processing Male

EXTENDED RANGE



Red
1.0 Lbs



Orange
2.0 Lbs



Green
4.0 Lbs



Gray
0.0 Lbs

STANDARD



Blue
1.5 Lbs



Pink
3.0 Lbs



Clear
5.0 Lbs



Purple
8.0 Lbs

Decreased Retention

Medium Retention

Increased Retention



8587



8586



8588



8589



8922
Locator Spot Face
Diamond Bur



8924 Pilot Drill

Locator Insert Driver (Latch 23mm)



Locator Insert Driver
(Square 23mm)




***See Abutment Listing for ALL Locator Abutments**



0700 0217
CM LOC® Extractor for
female part housing.
**compatible with Locator.*

- 8390 Locator Core Tool (End) Final
- 8393 Locator Core Tool
- 8394 Locator Abutment Holder Sleeve
- 8397 Locator Male Removal End
- 44-TW30 30 Ncm Torque Wrench
- 44-TW20 20 Ncm Torque Wrench



Abutment Driver

- 8505 Locator Impression Coping (4)
- 8505-20 Locator Impression Coping (20)
- 8516 Locator 5mm Ø Female Analog (4)
- 8516-20 Locator 5mm Ø Female Analog (20)
- The 5mm Ø analog is for Wide Platform Implants and Natural Teeth*
- The 4mm Ø analog is for Narrow and Regular Platform Implants*
- 8530 Locator 4mm Ø Female Analog (4)
- 8530-20 Locator 4mm Ø Female Analog (20)



8505 Impression Coping



8530 Female 4mm Analog

- 8517 Locator Parallel Post (4)
- 8514 Locator Block Out Spacer (20)
- 8569 Locator Processing Spacers (4)
- 8569-20 Locator Processing Spacers (20)
- 8519-2 Locator Male Processing Pkg TI (2)
- 8519-10 Locator Male Processing Pkg TI (10)
- TI housing, blockout spacer, blue / pink / clear males*



8519-2 Male Processing Package

- 8540-2 Locator Extended Range Male Proc Pkg (2)
- 8540-10 Locator Extended Range Male Proc Pkg (10)
- TI housing, blockout spacer, red / orange / green males*



8540-2 Ext Range Male Processing Pkg

- 8550-2 Locator Male Processing Pkg SS (2)
- 8550-10 Locator Male Processing Pkg SS (10)
- Stainless Steel housing ideal for welding to CoCr framework*
- SS housing, blockout spacer, blue / pink / clear males*



8028-2 Bar Processing Package

- 8028-2 Locator Bar Processing Package (2)
- 8028-10 Locator Bar Processing Package (10)
- TI hous. w/ yellow, 2 block. spacers, blue / pink / clear*



8520 8521 8522 8528

- 8520 Locator Root 0 Degree Female (2)
- 8520-10 Locator Root 0 Degree Female (10)
- 8521 Locator Root 10 Degree Female (2)
- 8521-10 Locator Root 10 Degree Female (10)
- 8522 Locator Root 20 Degree Female (2)
- 8522-10 Locator Root 20 Degree Female (10)
- 8528 Locator Root Cast-to Female (2) (Coping)
- 8528-10 Locator Root Cast-to Female (10) (Coping)



8519 - LOCATOR MALE PROCESSING PACKAGE TI

Includes: **Titanium** Denture Cap with Black Processing Male, White Nylon Male (5lbs), Pink Medium Retention Nylon Male (3lbs), Blue Light Retention Male (1.5lbs), and Block-Out Spacer.



8550 - LOCATOR MALE PROCESSING PACKAGE SS

Includes: **Stainless Steel** Denture Cap with Black Processing Male, White Nylon Male (5lbs), Pink Medium Retention Nylon Male (3lbs), Blue Light Retention Male (1.5lbs), and Block-Out Spacer. **Designed for laser welding to cast frameworks.**



8540 - LOCATOR EXTENDED RANGE MALE PROCESSING PACKAGE TI

Includes: **Titanium** Denture Cap with Black Processing Male, Red Male (0.5-1.0lbs), Orange Male (2.0lbs), Green Retention Male (3.5-4.0lbs), and Block-Out Spacer.



8505 - LOCATOR IMPRESSION COPING



8530 - LOCATOR FEMALE ANALOG 4MM NARROW AND REGULAR PLATFORM



8516 - LOCATOR FEMALE ANALOG 5MM WIDE PLATFORM