


Hader Bar System

45-531060 Hader Alignment Housing




6pcs

Alignment Housing



Incorrect - Female clip cannot flex over the bar

- 45-531010 Hader Yellow Riders- Standard
- 45-531012 Hader White Riders- Decreased
- 45-531014 Hader Red Riders- Increased
- 45-0101760 Hader Blue Riders- Worn Bars



- 45-0101775 Hader Aluminum Analog
- 45-0101703 Hader Green Processing Spacers

- 45-0101704 Hader Plastic Bar Patterns (Standard)
- 45-0101709 Hader Plastic Bar Pattern (Short)



Short & Standard Hader Bar Patterns



Delrin Analog Bar



Aluminum Analog



Hader Insertion Tool
45-0101705

Metal Riders for Hader/Dolder Bars/Round Bars

Metal Riders provide both rotational and vertical movement with a choice of two styles of acrylic retentions

- 40-05050010 Ackermann Clip A - BL
- 40-05050011 Ackermann Clip B - MD



- 40-050527 CM Rider Clips (single)
- 40-055801 CM Rider Clips (5 pack)

- 40-05000679 Round Bar with Rider E (50mm)

- 40-054747 Dolder Macro E (25mm)
- 40-052046 Dolder Macro E (50mm)
- 40-054746 Dolder Micro E (25mm)
- 40-052043 Dolder Micro E (50mm)



- 40-05000681 Dolder Macro T (47.5mm)
- 40-05000680 Dolder Micro T (47.5mm)

Universal Plunger Loc

Universal size Titanium Locking Plunger for Overdenture Bars

7.5mm long, may be cut to 6, 5, or 4mm. 1.5mm \varnothing plunger
Plunger Diameter compatible with Swiss-Loc & Lew Passive

- 1776000 Universal Plunger Loc Attachment
- 1776005 Universal Plunger Loc Processing Jig
- 1776010 Universal Plunger Loc Male Only
- 1776030 Universal Plunger Loc Delrin inserts(ea)
- 1776015 Universal Plunger Loc Twist Drill

