

# OT Equator Components



VIOLET CAP  
STRONG RETENTION  
Ref. 140CEV



WHITE CAP  
STANDARD RETENTION  
Ref. 140CET



PINK CAP  
SOFT RETENTION  
Ref. 140CER



YELLOW CAP  
EXTRA SOFT RETENTION  
Ref. 140CEG



BLACK CAP  
PROCESSING  
Ref. 140CEN



Ref. 041CAN B



Ref. 041CAM B

STAINLESS STEEL  
HOUSING



Ref. 055SCN



Ref. 055SCM

OT EQUATOR PROFILE  
CASTABLE  
MALE ATTACHMENT

**Components for OT STRATEGY  
also available**

# OT CAP System Components

available in normal 2.5mm and micro 1.8mm



CLEAR CAPS  
STANDARD RETENTION  
Ref. 040CRN - Ref. 040CRM



PINK CAPS  
SOFT RETENTION  
Ref. 040CRN SN - Ref. 040CRM SN



YELLOW CAPS  
EXTRA SOFT RETENTION  
Ref. 060CRN AY - Ref. 060CRM AY



GREEN CAP  
VERY ELASTIC RETENTION  
Ref. 049PCN - Ref. 049PCM



GREY CAP  
FIRM RETENTION  
Ref. 060CRN AC - Ref. 060CRM AC



LIGHT GREEN CAP  
STRONG RETENTION  
Ref. 040CRNDR8 - Ref. 040CRMDR8  
(undersized internal diameter for worn sphere)



ORANGE CAP  
VERY ELASTIC (light) RETENTION  
Ref. 049PCNDR8 - Ref. 049PCMDR8  
(undersized internal diameter for worn sphere)



BLACK CAP  
PROCESSING  
Ref. 043CLN - Ref. 043CLM  
(FOR LABORATORY USE ONLY)



## Cuff Height Measurer

The Cuff Height Measuring Gague from Rhein83 determines the exact tissue height requirement for any implant brand with an internal or external connection. Easy to read, the color-coded millimeter measurement tool accurately measures tissue height from 0.5 - 7mm.

*The kit contains one cuff height slider gauge, one cuff height fixed rod gauge, one silicon ring dispenser and 20 silicon rings.*



Ref. 008MBG



2 Transfer Pivots - Normal  
Ref. 044PPN



2 Transfer Pivots - Micro  
Ref. 044PPM

SPHERO Block Analog Kit  
Normal

Ref. 044ABN

SPHERO Block Analog Kit  
Micro

Ref. 044ABM