



# Sterngold™ ERA

The ERA attachment system is Sterngold's response to the dental professionals need for simple, dependable and highly versatile attachment system for removable prostheses. Sterngold is used in removable partial dentures treatment, tooth supported overdentures and implant cases, ERA is an excellent choice for both experienced practitioners and those new to attachment procedures.



## Stern ERA Micro Complete Attachment - Starter kit - 811012

### Extracoronary semi-precision attachment - 811011

Kit contents: 1 female, 2 black males, 2 white males, and 1 orange male

Micro ERA Starter Kit	811011
Micro Partial ERA Metal Jacket	811021
Micro ERA Partial Denture Impression Coping	811235
Micro Partial ERA Male Black - (5)	811013
Micro Partial ERA Male White (5)	811014
Micro Partial ERA Male Orange (5)	811015
Micro Partial ERA Male Blue (5)	811016
Micro Partial ERA Male Grey (5)	811017
Micro Partial ERA Male Yellow (5)	811018
Micro Partial ERA Male Red (5)	811019
Micro Partial ERA Male Assorted 7	811009
Micro ERA Female	811010
Micro ERA Partial Denture Processing Jig	811020
Micro ERA Laboratory Tool Kit	811025
Micro ERA Dentist Tool Kit	811026
Micro ERA Core Cutter Bur	811023
Micro ERA Seating Tool	811022
Micro ERA Paralleling Mandrel	811024
EZ Pickup™ Syringe and Tips (Self-Curing)	220235
ERA Attachment Extraction Tool	811027

## Stern ERA-RV Partial Attachment - Starter kit - 811115

### Resilient semi-precision attachment

Contents: 1 female and 2-Black males 2-white males, and 1 orange male

Stern ERA-RV Starter kit	811115
ERA-RV Metal Jacket - Partial	811168
ERA RV Impression Coping	811232
ERA-RV DE Female SS	811610
ERA Male Black RV Partial 5	811125
ERA Male White RV Partial (5)	811135
ERA Male Orange RV Partial (5)	811145
ERA Male Blue RV Partial (5)	811155
ERA Male Grey RV Partial (5)	811165
ERA Male Yellow RV Partial (5)	811175
ERA Male Red RV Partial (5)	811176
ERA Male Assorted RV Partial (7)	811166
ERA Female only	811190
ERA/ERA-RV partial denture processing jig	811250
ERA Laboratory Tool kit	811200
ERA Dentist Tool Kit	811240
ERA core cutter bur	811220
ERA seating tool	811230
ERA paralleling mandrel	811210
EZ Pickup™ Syringe and Tips (Self-Curing)	220235
ERA Attachment Extraction Tool	811027



ERA-RV  
Metal Jacket - Partial  
811168



ERA RV  
Impression Coping  
811232


ERA/ERA-RV  
partial denture  
processing jig  
811250

## Stern ERA Overdenture Attachment - 811310

### Resilient semi-precision attachment

Contents: 1 female and 2-Black males 2-white males, and 1 orange male

Micro ERA Overdenture Starter Kit	811045
Micro ERA Metal Jacket Overdenture	811043
Micro ERA Overdenture Impression coping	811236
Micro ERA Male Black Overdenture, 5	811035
Micro ERA Male White Overdenture 5	811036
Micro ERA Male Orange Overdenture 5	811037
Micro ERA Male Blue Overdenture 5	811038
Micro ERA Male Grey Overdenture 5	811039
Micro ERA Male Yellow Overdenture 5	811040
Micro ERA Male Red Overdenture 5	811041
Micro ERA Male Assorted Overdenture 7	811029
Micro ERA Overdenture Female	811030
Micro ERA Overdenture Processing Jig	811042
Micro ERA Laboratory Tool Kit	811025
Micro ERA Dentist Tool Kit	811026
Micro ERA Core Cutter Bur	811023
Micro ERA Seating Tool	811022
Micro ERA Paralleling Mandrel	811024
EZ PickUp™ Syringe and Tips (Self-Curing)	220235
ERA Attachment Extraction Tool	811027



The six different males provide consistent retention throughout the life of the attachment. ERA Males use the same colour code.

#### Black is the processing male

White = smallest button diameter (least retention)

Orange

Blue

Grey

Yellow

Red = largest button diameter (most retention)



## Stern ERA Overdenture Attachment - 811310

### Resilient semi-precision attachment

Contents: 1 female and 2-Black males 2-white males, and 1 orange male



#### ERA PickUp Syringe and Tips (63-220235)

1 Syringe @ 15 grams and 15 tips

#### SternVantage Varnish LC (63-221001) 30 ml bottle

ERA PickUp will chemically bond to acrylic only if a primer is used. Use SternVantage Varnish LC as a primer for acrylic.

Stern ERA Overdenture Starter Kit	811300
ERA Metal Jacket Overdenture	811380
ERA Overdenture Impression Coping	811233
ERA Male Black Overdenture (5)	811320
ERA Male White Overdenture (5)	811330
ERA Male Orange Overdenture (5)	811340
ERA Male Blue Overdenture(5)	811350
ERA Male Grey Overdenture (5)	811360
ERA Male Yellow Overdenture (5)	811370
ERA Male Red Overdenture (5)	811375
ERA Male Assorted Overdenture (7)	811365
ERA Overdenture female	811390
ERA Overdenture processing jig	811395
ERA Laboratory Tool kit	811200
ERA Dentist Tool Kit	811240
ERA core cutter bur	811220
ERA seating tool	811230
ERA paralleling mandrel	811210
EZ PickUp™ Syringe and Tips (Self-Curing)	220235
ERA Attachment Extraction Tool	811027